# ventura flow

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## EN INSTRUCTIONS FOR USE

### 1. DEFINITION

**Ventura Flow** is a light cure easy flowing microhybrid composite for anterior and posterior restorations, and for sealing. It is based on methacrylate-resin and inorganic filler particles of 0,05-2,6 µm. The total filler load is 53,00 %, the total filler volume is 28,00 %. The composite material meets the EN ISO 4049 requirements.

### 2. COMPOSITION

**Ventura Flow** contains:

Resina based on BIS-GMA, radio-opaque glass filler silanized, catalysts, inhibitors, pigments.

### 3. INDICATIONS

A. Class I, II, III, IV and V filling according to Black.

B. Sealings.

C. Reconstructions of acrylic and ceramic pieces.

### 4. COSMETICAL AND PHYSICAL PROPERTIES

Very aesthetic.

Very high strength.

Abrasion resistance.

Strong radiopaque.

Excellent colour stability.

High gloss polishability.

Highly filled.

Low volumetric shrinkage.

Extremely low contraction stress.

### 5. CONTRAINDICATIONS

Exposed pulp, pulpitis, known allergy to methacrylate resins.

### 6. SIDE EFFECTS

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore, for pulp protection, the cavity floor, in cases of deep excavations, should be covered with a thin layer of hard setting calcium hydroxide material.

**Ventura Flow** does not contain any TEGDMA or HEMA

### 7. INCOMPATIBILITY WITH OTHER MATERIALS

Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the composite and causes discoloration. Do neither store the composite material in proximity of eugenol containing products, nor let the composite allow to come into contact with materials containing eugenol.

### 8. PREPARATION OF TOOTH AND COLOUR MATCHING

Clean the tooth with a fluoride-free polishing paste prior to preparation and colour matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate according to VITA® -Shades.

### 9. APPLICATION

After isolation, cavity preparation, pulp protection (in case of deep excavations), etching and application of a light cure adhesive system,

**Ventura Flow** can be applied.

### 10. PLACEMENT OF VENTURA FLOW

Apply **Ventura Flow** in the selected shade. Place it by using a non-metallic instrument. Transparent matrix strips may be used.

### 11. CURING OF VENTURA FLOW

Place increments of no more than 2 mm. Light cure each layer for 40 seconds.

The indicated light curing time refers to halogen- or LED-lamps with 500-600 mW/cm<sup>2</sup> and 400-500 nm. It is valid for an increment thickness of maximum 2 mm and for all shades. When using matrices, after removing them expose filling to light again from lingual and buccal.

light power	≥ 500 mW/cm <sup>2</sup>	≥ 800 mW/cm
colour opaque	40 sec.	30 sec.
all other colours	40 sec.	20 sec.

### 12. FINISHING

Trim filling with 40 µm and 12 µm diamonds, points or burs. Polish to high gloss with discs or silicone points. Interproximal finishing is accomplished by fine grit finishing strips.

### 13. MARKING OF THE PACKAGING

The lot number and expiry date are printed on each package.

### 14. STORAGE

Do not store above 25 °C (78 °F). Prolonged temperature above 28 °C can shorten the shelf life of the product.

### 15. EMERGENCY PROCEDURES

Contact with the skin : wash with soap and water.

Contact with the oral mucous membrane : rinse with large quantities of water.

Contact with the eyes : wash immediately with large quantities of water and consult an ophthalmologist.

### 16. HYGIENE

Use the application instruments for only one patient. Clean and disinfect all other instruments after each use. To avoid contamination, dispense the materials at a distance from the patient.

### 17. ADDITIONAL NOTES

Do not use any resin to adjust viscosity of composite restorative material. Contact of resin pastes with skin should be avoided, especially by anyone having known resin allergies.

VITA is a registered trade mark of the VITA-Zhanfabrik, Bad Säckingen, Germany.

Instructions for use issued in October 2016 (2016-10).

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