# ventura bite

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## EN INSTRUCTIONS FOR USE

**Ventura BITE** is a silicone elastomer of the vinyl-polysiloxane type whose reticulation takes place based on a polyaddition reaction. The reticulated product's mechanical properties and reproduction of detail are excellent.

**Ventura BITE** presents a high tixotropy and low contraction. It does not adhere to organic resins or gypsum.

### INDICATIONS

Bite Registration

### PREPARATION OF THE DISPENSING GUN

1. Load the cartridge by sliding it vertically through the leads provided for this purpose.
2. Press the lever until the plunger rests on the cartridge.

### MIXING

Unscrew the cartridge cap by turning it 90° (1/4 of a revolution), remove it and discard it, as it cannot be reused. Extrude impression material until it exits both chambers at the same rate and then mount the appropriate mixing tip, blocking it by means of a 90° turn (1/4 of a revolution).

### MAKING THE IMPRESSION

Extrude material first onto areas where teeth have been prepared and then onto the rest to areas and ask the patient to bite. The material should be placed directly on the teeth within 30 seconds from starting the mix and kept in position for at least 45 seconds. Once you have finished the procedure, leave the mixing tip on the cartridge until the next use as it serves as a cap. Between uses, store the dispensing gun with the mixing tip facing down. For the next mix using the same cartridge, remove the used mixing tip by turning it 90° (1/4 of a revolution), insert a new mixing tip and proceed.

### DISINFECTING

Once bite registration is removed from mouth, rinse under cold water. The impression can then be disinfected with appropriate disinfecting solutions or sprays.

### PRECAUTIONS

In sensitive subjects, vinyl polysiloxanes may cause irritation or other allergic reactions.

Contact with latex gloves may prevent vinylpolysiloxane from setting. Avoid direct or indirect contact with latex gloves. Instead, wash your hands thoroughly to eliminate all the traces of impurities or use vinyl gloves.

Silicone adheres to clothes from where it may be very difficult to remove, even by using chemical agents.

### STORAGE

In order for the product to maintain its optimum characteristics, it is advisable to store it at temperatures between 10°C and 25 °C (50°F-77°F) and use it before the expiry date.

### LOT NUMBER AND EXPIRY DATE

They are printed on each package.

### TECHNICAL CHARACTERISTICS

Classification ISO 4823	: Type 2 - Medium bodied consistency.
Color	: Blue
Mixing ratio base : catalyst	: 1 : 1
Mixing time	: Automixing
Time in oral cavity	: 30" *
Working time at 23 °C	: 45" *
Setting time	: 1'15" *
Linear dimensional variation at 24 hours	: < 0.02 %
Shore A Hardness at 24 hours	: 93.

\*These times are intended from the beginning of mixing at 23°C (73°F). Higher temperatures shorten, lower temperatures prolong these setting times.

Instructions for use issued in September 2013 (2013-09)

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