# ventura dual core built

automix

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## EN INSTRUCTIONS FOR USE

**Dual curing flowable hybrid composite material for core build-ups and post cementation.**

The auto-mixing **Ventura dual core built** is a radio-opaque dual curing hybrid compound for core build-ups and post cementation. It is available in self-mixing MINIMIX 1:1 syringes, which allow direct application. The product has a short setting time without producing a large amount of heat. Its dual cure feature allows the dentist to carry out core build-ups and post cementations also in those cases in which it cannot be guaranteed that light-curing would be sufficient. **Ventura dual core built** is composed of methacrylate resins and 0,05-1,5 µm inorganic filler particles.

The filler quantity is 66,00% by weight and 46,00% by volume.

**Ventura dual core built** complies with the requirements of the norm EN ISO 4049 standard.

### 1. Preparing the MINIMIX syringe

**Ventura dual core built** is mixed automatically as it is dispensed applying a smooth uniform pressure.

Remove the MINIMIX syringe cap and throw it away (do not reuse it). Instead, place the special 1:1 mixing cannula. Turn the cannula 90° until it is locked in position. Dispose of the first 2-3 millimetres of material to ensure a homogenous mix. At this moment, the material is ready to be applied.

Note: Store the syringe with the used mixing cannula in a dark place.

### 2. Working and curing time

At room temperature (23°C/73,4°F)

Working time	1'30"
Self-polymerisation time	3'30"±10"
Photo-polymerisation time	40"

Light-cure each layer with a halogen light bulb (=350mW/cm²) or an LED (=850mW/cm²). Wavelength for correct polymerisation: 470±20nm.

### 3. Post cementation

#### 3.1. Insulation

A rubber dam is recommended to insulate the tooth.

#### 3.2. Preparing the radicular canal

Follow the instructions of the post manufacturer.

Eliminate all old restorations and decay of the tooth. Prepare and clean the root canal with, e.g., a sodium hypochlorite solution, and rinse and remove excess solution from the canal with a paper point.

#### 3.3. Applying an adhesive agent

Condition with orthophosphoric acid and apply a dual-cure adhesive agent following the manufacturer's instructions. Otherwise, use a self-etching adhesive agent.

#### 3.4. Post cementation

Prepare the selected post in accordance with the manufacturer's instructions.

Apply **Ventura dual core built** into the prepared root canal and to the post. Carefully place the post in the canal and keep a firm pressure for 5-10 seconds until the post is well seated. To stabilise the post, light-cure the crown part of the cemented post for 40 seconds, placing the light source as close as possible to the post. Self-curing takes place in 3:30 minutes.

Once **Ventura dual core built** has set, proceed with the reconstruction process.

### 4. Core Build-Up

#### 4.1 Insulation

A rubber dam is recommended to insulate the tooth.

#### 4.2. Preparing the cavity

Eliminate all old restorations and decay of the tooth. If necessary, place pins or posts. Follow the post manufacturer's instructions.

#### 4.3. Pulp Protection

Should it be a vital tooth, when deep excavations are going to be made, cover the cavity floor with a thin layer of quick setting calcium hydroxide.

#### 4.4. Applying an adhesive agent

Condition with orthophosphoric acid and apply a dual-curing adhesive agent following the manufacturer's instructions. Otherwise, use a self-etching adhesive agent.

Note: It is essential that the etched surfaces of dentine and enamel are dry and free from contamination before applying **Ventura Dual Core Built**.

#### 4.5. Applying the composite

Place the mixing cannule directly in the prepared part and press to extract the mixture. The cavity is filled bottom upwards to avoid leaving voids. To simplify the application of the product, place a dental matrix band around the prepared tooth.

**Ventura dual core built** can be modelled using instruments for composites. An flat-ended interproximal carver is recommended.

Apply **Ventura dual core built** directly into the prepared tooth and allow the system to self cure for 3:30 minutes. After that time, the material must be light-cured (40 seconds), placing the light source as close as possible to the restoration. This technique will ensure that the best physical properties for the material are obtained.

A dental probe can be used to check whether **Ventura dual core built** has set completely. Do not remove the matrix until the material has set.

Final crown preparation can be carried out using crown preparation burs.

### 5. Further note

Do not use resins to adjust the viscosity of the reconstruction material.

Do not use this product combined with substances that may contain eugenol, since eugenol may affect the setting of the composite and cause discolouration. Store the material away from products containing eugenol and avoid the material coming into contact with products containing eugenol.

Contact between the resin pastes and the skin and gingival tissue must be avoided, especially when there is a known allergy to resins.

### 6. Storage

Do not store at temperatures above 12 °C. Keep unopened material in the refrigerator. Opened cartridges must be used within 3 months

Do not use the product after its expiry date.

### 7. Contraindications

Allergy to methacrylates. Open pulp or pulpitis

Keep away from children

Only for dental use

Instructions for use published in 2017-12 (December 2017).

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